

## ABSTRACT FOR 2010 NSMMS

### Improved Creep Measurements for Ultra-High Temperature Materials

Robert W. Hyers<sup>1</sup>, X. Ye<sup>1</sup> and Jan R. Rogers<sup>2</sup>

<sup>1</sup> University of Massachusetts, Amherst

<sup>2</sup> NASA/MSFC/EM50

Our team has developed a novel approach to measuring creep at extremely high temperatures using electrostatic levitation (ESL). This method has been demonstrated on niobium up to 2300°C, while ESL has melted tungsten (3400°C). The method has been extended to lower temperatures and higher stresses and applied to new materials, including a niobium-based superalloy, MASC.

High-precision machined spheres of the sample are levitated in the NASA MSFC ESL, a national user facility, and heated with a laser. The samples are rotated with an induction motor at up to 30,000 revolutions per second. The rapid rotation loads the sample through centripetal acceleration, producing a shear stress of about 60 MPa at the center, causing the sample to deform. The deformation of the sample is captured on high-speed video, which is analyzed by machine-vision software from the University of Massachusetts. The deformations are compared to finite element models to determine the constitutive constants in the creep relation. Furthermore, the non-contact method exploits stress gradients within the sample to determine the stress exponent in a single test.





# Improved Creep Measurements for Ultra-High Temperature Materials

Robert W. Hyers<sup>1,2</sup>, X. Ye<sup>1</sup>, and Jan R. Rogers<sup>3</sup>

<sup>1</sup>University of Massachusetts Amherst, <sup>2</sup>R. Hyers and Associates, Amherst, MA

<sup>3</sup>NASA Marshall Space Flight Center

## Motivation


- Increasing need for high-temperature materials
- Higher operating temperatures leads to greater performance and efficiency
- Creep of metals is important at high temperatures ( $T \geq 0.4-0.5 T_{\text{melt}}$ )
- High-temperature materials ( $T_{\text{melt}} \geq 2500^\circ\text{C}$ ) are being developed and ready to use
  - i.e. ultra-high temperature ceramics and platinum group metals
- Conventional methods limited to  $\sim 1700^\circ\text{C}$
- Non-contact method demonstrated up to  $2350^\circ\text{C}$

## Applications

- Next Generation turbine blades
  - $>1250^\circ\text{C}$  for more than 4000 hours
- Rocket Nozzle
  - Up to  $3000^\circ\text{C}$ , high stress
- Hypersonic Flight
  - Leading edge materials
  - $\geq 2700^\circ\text{C}$


## Conventional Creep Tests

- Specimen in contact with test equipment
- Materials become reactive at high temperatures and incompatible with the containers or equipment




## Non-Contact Creep Tests

- 2-3 mm diameter high-precision spheres
- Load by centripetal acceleration,
- Rotation rates up to 30,000 rev/sec
- Loads up to 100 MPa, Temperature to  $2300^\circ\text{C}$




## Non-Contact Creep Tests



## Non-Contact Creep Test Analysis

- Deformed shape depends on stress exponent
  - The ratio of the polar to equatorial strains (Strain Ratio) from a single ESL test determines the stress exponent
- FEA model used to generate a stress exponent versus strain ratio plot




## Induction Motor Design

- Clearance: unobstructed view of levitated sample.
- Vacuum compatibility: Materials, cooling.
- Integration with MSFC ESL
- Performance:
  - up to 30,000 rev/sec, 10x better
  - up to  $\sim 5 \times 10^{-10}$  N-m, 10x better
- Lower temperature materials
- Higher speed = higher stress
- Sponsored by NASA IPP
- Completed July 2009

## Rotation by Induction Motor

Decouple temperature and stress

- Increased acceleration  $\Rightarrow$  greater maximum stress
- Enables experiments below  $2000^\circ\text{C}$
- Torque depends on electrical conductivity



## Current Work

- Non-contact measurement of creep in ZrB<sub>2</sub> and ZrB<sub>2</sub>-SiC composites.
- Material from NSWCCD
  - ZrB<sub>2</sub>: 6 mm grains
  - SiC: 2 mm grains
- Spheres machined ITI
- Tests at NASA MSFC
- Analysis with UMass
- Comparison to flexural creep measurements





Fig. 7. Stress exponent for creep deformation of ZrB<sub>2</sub>/SiC ceramics at 1400°C as a function of SiC content and particle size (2 and 10 μm).

## Status: Modeling and Analysis

- FE Model running with parameters extrapolated from Talmy, *et al.*
- Pure ZrB<sub>2</sub>: Model predicts 100 hours to 10% strain at  $2000^\circ\text{C}$  and 100 MPa (rotation rate 32,500 Hz). 2.8 GPa / 150,000 Hz needed for 2 hour experiment.


- ZrB<sub>2</sub> + 25 vol% SiC: Model predicts 1.5 hours to get 10% equatorial strain at  $1900^\circ\text{C}$  and 100 MPa (rotation rate 32,500 Hz).



X. Ye and R.W. Hyers, "Computational Simulation of Non-contact Creep Deformation of ZrB<sub>2</sub>/ZrB<sub>2</sub>+SiC", *Journal of the European Ceramic Society* 30, 2191-2196 (2010).

## ZrB<sub>2</sub> and ZrB<sub>2</sub> + 25 vol% SiC

- ZrB<sub>2</sub> + 25 vol% SiC at  $1800^\circ\text{C}$
- Samples accelerated by ESL induction motor
  - 22,200 Hz (1.3 million RPM): 232 minutes at load, temperature
  - 25,400 Hz (1.5 million RPM): 219 minutes at load, temperature
- Total strain  $\sim 2\%$ , 2X less than extrapolation
- Pure ZrB<sub>2</sub> measured at  $2000^\circ\text{C}$ 
  - 27,100 Hz (1.6 million RPM), 333 minutes at load, temperature.
- Total strain  $< 0.07\%$ : consistent with extrapolation



## On-Going Work

- Continue work with UHTC's, Ni- and Nb-based superalloys, other materials.
- Improved treatment of evaporation
- Further analysis of non-contact creep method
  - Multiple creep mechanisms
- Higher stresses
  - Photon pressure:  $\sim 3,000$  rev/sec  $\Rightarrow \sim 1$  MPa max.
  - Present induction motor and measurement:  $\sim 30,000$  rev/sec  $\Rightarrow \sim 100$  MPa max.
  - NASA developing new rotation measurement: even higher speeds, stresses.
- Higher torque to reduce experiment duration.
- More automation of analysis.

2010 NSMMS

This work was supported by the NASA Innovative Partnerships Program and AFOSR STTR Grant FA9550-09-C-0089